

Please contact Invacare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
Or visit our website : www.invacare.eu.com

Assembly by trained
service providers or
Invacare



alber® Viamobil® V25

Adaptable electric motor pushing assistance device with attendant control only



Description:

- Full motorization for pushing and braking with ergonomic control handle
- Wheelchair charger and attachment kit
- 6.6 Ah battery.
- 2-year warranty including battery excluding wearing parts.
- Maximum user weight: 160 kg



Date :
Purchase Order :
Customer Code :

Distributor Stamp:

Indicative Public prices - Use the online configurator on Invacare Pro to generate quotes or orders free of incompatibilities with the last rates in force

Please make your selection with a cross in the box ◆ Standard feature ○ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros
Viamobil®V25				
1521003	<input type="checkbox"/>	alber Viamobil® V25		4 197,16
MANDATORY OPTIONS				
1457771	<input type="checkbox"/>	Bracket for mounting to the wheelchair		236,97
1584856-VI	<input type="checkbox"/>	Installation of brackets by Invacare (return transport costs not included). Average assembly time: 5 working days upon receipt of the chair (In the case of an existing chair, plan to return it to our workshops within 10 days).		214,22
RECOMMENDED OPTION				
1584856-2VI	<input type="checkbox"/>	Installation of bars Viamobil V25 (when the brackets are selected in mounted option on the order form of Action et Rea wheelchairs)		189,57

To determine compatibility:

Invacare® wheelchair (Action®, Rea®, kuschall®), please indicate the serial number of the wheelchair to be equipped:

Wheelchair of all other brands, please send us with this sheet the measurement of the armchair to be equipped.

Please note: Wheelbase of the chair subject to change. Avoid all options close to the rear bearings which could interfere with assembly.

OPTIONS

1593205	<input type="checkbox"/>	Anti-tipping castors, usable only when the motorization is removed from the wheelchair		355,45
1593418	<input type="checkbox"/>	Removable additional handles		229,38
1520929	<input type="checkbox"/>	Tilt tube (only usable when the motorization is removed from the chair)		113,74
1593389	<input type="checkbox"/>	Viamobil® anti-tip wheels, pair, mandatory for slopes >8% (to be installed on the motorisation)		355,45



1593389



1593205



1593418



1520929

alber® Viamobil® V25 Collectivité

Adaptable electric motor pushing assistance device with attendant control only



Description:

- Pushing and braking assistance motorization with height-adjustable handles and anti-tip wheels as standard.
- To be combined with universal mounting for Rea Clematis Pro & Clematis Pro E-tilt
- Requires the removal of the chair's original anti-tippers
- NEW: Allows the use of a single motorization for Clematis Pro & Clematis Pro E-tilt with different seat width and seat height.
- 6.6 Ah battery / 2-year warranty included battery excluding wearing parts.



Date :
Purchase Order :
Customer Code :

Distributor Stamp:

Please make your selection with a cross in the box ◆ Standard feature ○ Free of charge option

CODE	Qty	DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros
Assembly of the viamobil® set				
1521003 + 1593418 + 1593389	<input type="checkbox"/>	alber Viamobil®V25 Collectivité (V25 motor 1521003, additional removable handles 1593418, Viamobil V25 anti-tip wheels 1593389 and fixing bars included)		4 781,99
MANDATORY OPTIONS				
1593893 (width 390) 1594176 (width 440) 1594177 (width 490)	<input type="checkbox"/>	Clematis Pro & Clematis Pro E-tilt universal brackets, unit, provide 1 per chair NB: - Increases overall width by 20mm on Clematis Pro seat width 390mm - Requires the removal of the original anti-tippers		418,01
1584856-VI	<input type="checkbox"/>	Installation of brackets by Invacare, unit, provide 1 per chair (return transport costs not included) Average assembly time: 5 working days upon receipt of the chair (In the case of an existing chair, plan to return it to our workshops within 10 days).		214,22
OPTIONS				
1593205	<input type="checkbox"/>	Anti-tipping castors, usable only when the motorization is removed from the wheelchair		355,45